



AdultLife Ventilator

saving lives is our mission

Introduction



- The AdultLife Ventilator is a pressure controlled, time cycled ventilator and CPAP device. Its simple yet effective and durable design allows for affordability, ease of use, and portability.
- AdultLife's Core Features
 - Internal Battery
 - Internal Air Pump
 - Low/High/Disconnect Pressure Alarm
 - Low Maintenance
 - Low Cost



Specifications



- Peak Inspiratory Pressure (PIP)
 - 0 to 40cm H₂O
- Positive End Expiratory Pressure (PEEP)
 - 0 to 20cm H₂O
- Inspiratory Time (IT)
 - 0 to 2.5 Seconds
- Expiratory Time (ET)
 - 0 to 5 Seconds
- Breath Rate Per Minute (BPM)
 - 0 to 35 Breaths/Min
- Tidal Volume
 - 200-1000cc
- Oxygen Input/Medical Air Inputs
- Optional: Internal Pump (60L/min)
- Optional: Internal Battery (2-12 hr)
- Compatible with 22mm corrugated tubing and accessories (filter, humidifier, etc.)
- Internal Pressure Release Valves
 - Set at 40cm H₂O
- Alarms:
 - Disconnect/Power Off Alarm
 - Max PIP Alarm
 - PEEP Alarm
 - Occlusion Alarm
 - Hypoventilation
- 110-240V Compatible (US Type B plug)



FDA Compliance



- 3/23/2020 Proposed Emergency Use Approval (EUA) to FDA
- 3/25/2020 FDA Reviewed/Approved EUA Proposal and Requested Full Submission
- 4/6/2020 FDA EUA Submission*
- 6/17/2020 FDA EUA Approval**
- 10/15/2020 Estimated 510K Submission



* See appendix A for list of testing standards and medical review team

** <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/ventilators-and-ventilator-accessories-euas>

Training Support



- The AdultLife Ventilator includes the following resources for training and support:
 - Detailed instructions for use (IFU) available in both English and Spanish
 - Available training videos
 - USA respiratory university professor instruction on how to properly operate AdultLife ventilators
 - Availability to schedule video chat for further instruction with AdultLife medical staff and/or engineers



AdultLife Applications



- Respiratory Failure
- Respiratory Insufficiency
- Pneumonia
- Acute Respiratory Distress Syndrome (ARDS)
- Asthma
- Sepsis
- Brain Injuries
- CNS/Neurologic Conditions
- Trauma
- Drug Overdose



Warranty

- AdultLife ventilation devices are provided with a 1 year or total duration of the FDA defined Emergency Use Authorization (EUA) period, whichever comes first. If the EUA is terminated by the FDA, the 1 year warranty will continue upon FDA 510k approval. The components below are covered on the following schedule:
 - Internal Pump: 1 year or 10,000 hours
 - Pressure Gauge: 1 year
 - Internal Batteries: 1 year or 300 cycles



AdultLife Options



Base Model

AdultLife Ventilator
Wall Air/Cylinder Connectors
Patient Tubing
User Manual
Training Resources
1 Year Warranty

Complete

AdultLife Ventilator
Internal Pump
Internal Battery
Wall Air/Cylinder Connectors
Patient Tubing
User Manual
Training Resources
1 Year Warranty

Complete Plus

AdultLife Ventilator
Digital Monitoring System
Internal Pump
Internal Battery
Wall Air/Cylinder Connectors
Patient Tubing
User Manual
Training Resources
1 Year Warranty



Digital Monitoring System




A digital monitoring system is available with the AdultLife Pro Ventilator - Complete Plus.

This monitoring system allows for digital reading of tidal volume, pressure range, and wave form data. The system includes the following components:


- Digital Monitor
- Computer + Software Package
- Pressure Sensor System




Lung Model
Piston



Lung Model
Piston



Airway
Muscle
Trach./Alv.



Manufacturing Capabilities



- Manufactured in USA in compliance with ISO 13485 and 21 CFR part 802
- All components are sourced in the USA
- Limited quantities are available in stock for immediate fulfillment, please contact your representative for details.
- Orders under 500pcs can be fulfilled within 45 days.
- Orders under 1000pcs can be fulfilled within 60 days.





**AdultLife
Ventilator**

saving lives is our mission

Appendix A: Testing Standards and Medical Review

21 CFR Part 820 or ISO 13485: *Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes*

ISO 18562-1 First Edition 2017-03: *Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 1: Evaluation and Testing Within a Risk Management Process*

ISO 18562-2 First Edition 2017-03: *Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 2: Tests for Emissions of Particulate Matter*

ISO 18562-3 First Edition 2017: *Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 3: Tests for Emissions of Volatile Organic Compounds*

ISO 18562-4 First Edition 2017-03: *Biocompatibility Evaluation of Breathing Gas Pathways in Healthcare Applications - Part 4: Tests for Leachables in Condensate*

ISO 80601-2-12 First Edition 2011-04-15: *Medical Electrical Equipment - Part 2-12: Particular Requirements for the Safety of Lung Ventilators - Critical Care Ventilators [Including: Technical Corrigendum 1 (2011)]*

ISO 80601-2-80 First Edition 2018-07: *Medical Electrical Equipment - Part 2-80: Particular Requirements for Basic Safety and Essential Performance of Ventilatory Support Equipment for Ventilatory Insufficiency*

Additionally, we have had our medical team review for useability testing--something that is required for FDA 510K approval.

Medical Review Team:

Janelle Gardiner (Professor of Respiratory Care/Respiratory Therapist 15yrs exp.)

Phil Thaut (Lead Respiratory Therapist ICU--35+ yrs exp)

Clark Bishop (Pulmonologist)

Robert Guenter (Respiratory Therapist--former director of respiratory therapy at Utah Valley Hospital, 40+ yrs exp.)

David Eitel (Respiratory Therapist, current director of respiratory at Utah Valley Hospital, 15+ yrs exp.)