

Pumani *bubbleCPAP*

To order please contact us: info@hadleighhealth.com



The Pumani bubbleCPAP is used to treat respiratory distress syndrome (RDS) in premature neonates and respiratory illness in young children from 1 to 10 kg. Designed specifically for the low-resource hospitals, it is an effective, affordable and user-friendly solution currently in use in over 35 countries around the world.

Effective and Affordable - Provides the same level of therapeutic pressure as other commercial CPAPs for less than one-third of the price

User Friendly - Set-up training takes less than 5 minutes

Easy to Maintain - Minimal maintenance, spare parts kit included inside

Durable - Built and proven to run for years at a time

In Use In Over 35 Countries!

Pumani bubbleCPAP

To order please contact us: info@hadleighhealth.com

Oxygen Blending	21-90%
Gas Supply	Integrated Air Compressor
Total (blended) Flow Rate	0-10 L/min
Oxygen Flow Rate	0-10 L/min
Pressure	0-10 cm H2O
Voltage	220V and 50Hz
Power Consumption	15 W
Shipping Dimensions	48 cm x 38 cm x 33 cm
Shipping Weight	12kg
Regulatory Compliance	CE Mark
Accessories	Pumani Patient Circuit Kit, Bottles, Filters

