

Beluga 
resuscitator

LifeKit by 




BELUGA

An easy-to-use t-piece resuscitator suitable for emergencies and delivery units



Overview

- Portable
- Build-in battery
- Integrated blender
- Oxygen monitoring
- Easy-to-use touch screen
- Controlled target PIP and consistent PEEP
- No time restriction on use due to the gas tank limitation
- Usable in horizontal or vertical position



Physical Specifications

Unit dimensions	20 cm x 22 cm x 10 cm
Shipping (boxed) dimensions	30 cm x 25 cm x 15 cm
Net weight	1.5 kg
Shipping weight	2.0 kg

Performance Specifications

Power characteristics	100 - 240 VAC, 50-60 Hz
Manometer accuracy	± 0.5 cm H ₂ O
Peak Inspiratory Pressure (PIP)	5 - 50 cm H ₂ O
Positive End Expiratory Pressure (PEEP)	1 - 20 cm H ₂ O
Delivered oxygen concentration	21 - 100%
Integrated air compressor	Blower
Battery	Li-ion, 12V-4000mAh

Environmental Specifications

Operating	Ambient temperature: 0°C to +40°C Humidity: 0% to 90% RH non-condensing
Transport and Storage	Ambient temperature: -20°C to +60°C Humidity: 0% to 90% RH non-condensing

Standard Compliance

EN 60601-1:2006/A1:2013
EN 60601-1-2:2015

Ordering Information and Servicing

MAIN UNIT	ITEM ORDER NO
MTTS Beluga Resuscitator	RES-10000
Reusable tube	RES-10011
Single use PEEP valve with T-piece	RES-10012
Single use T-piece circuit*	RES-10020
Silicon mask - size S (42mm)	RES-10030
Silicon mask - size M (50mm)	RES-10031
Oxygen Hose	RES-10040
Arm Bag	RES-10050
Power cable - 12V Car Lighter*	COR-10006
*Optional	
WARRANTY	
Standard Warranty	1 Year Subject to MTTs policy

MTTS

MEDICAL TECHNOLOGY TRANSFER AND SERVICES Co., LTD
ISO 13485: 2016 Certified Company



VIA GLOBAL HEALTH
viaglobalhealth.com

Distributor: